

FUJIFILM Sonosite, Inc.
21919 30th Dr. SE
Bothell, WA
98021-3904
USA

Tel: 1.425.951.1200
Fax: 1.425.951.1201
www.sonosite.com

System Portal – Setup and Enablement Quick Guide

1. Remote software updates and system monitoring

Administrators can configure and connect the system to a cloud service that allows remote software updates and system monitoring. You can download updates from the cloud and install them directly on the system or upload system logs to the cloud that allow remote service personnel to troubleshoot issues. An internet connection is required, and your network firewall must be configured to **allow HTTPS connection on TCP port 443** to the following URLs:

- <https://system-portal.sonosite.com>
- <http://crl.microsoft.com>
- <https://system-portal-nosni.sonosite.com>
- <https://ffss.sac.keyscaler.io/>
- <http://www.microsoft.com/pkiops/crl>
- <https://sonositesystemportal.blob.core.windows.net>

Once the system is configured and connected to the cloud, the system notifies you of any available software updates or connection errors. The notifications appear in various places including the system menu, next to **Administration** and **System Information** in the system settings list, and on the administration settings page. You can also enable non-administrative users to download and install updates.

For more information, please review the 'Remote software updates and system monitoring' section in your device's User Guide.

2. Configuring the System Portal connection

1. Using your administrative login information, log into the administrative settings page.
2. Tap **System Portal Connection**.

FUJIFILM Sonosite, Inc.

21919 30th Dr. SE
Bothell, WA
98021-3904
USA

Tel: 1.425.951.1200

Fax: 1.425.951.1201

www.sonosite.com

3. On the next page, tap **Configure**.
4. On the configuration page, select the check box to enable connection to the feature.
5. If the system is not yet connected to a network, follow one of the procedures outlined in the **Connecting to the Network** shown in the device's User Guide.
6. To use a proxy connection, select the check box and enter the proxy server information in the provided fields. Server address and port number are required. This information is maintained if the proxy connection is turned off.
7. Tap **Save** to save the proxy settings and verify the connection. Tapping **Reset** will reset all the fields to the last time they were saved.

3. Downloading and Installing software updates

Administrative users can also perform the following procedure on the System Portal page, which is available from the Administration page. The system needs to be connected to the Internet, System Portal, and to AC power.

1. On the system information page, tap **Check Update**. The system checks for an available software update.
2. Tap **Download**. The system begins downloading the update. You can pause or cancel the download or continue using the system during the download process.
3. Monitor the download by checking the status on the page.
4. When the download is complete, verify that the system is ready for an update by making sure it is connected to a power supply and ending any exams.
5. Tap **Install**.

4. Uploading logs

To upload logs directly from the device, the system will need to be connected to the System Portal. Once the connection has been established, logs can be uploaded by performing the steps below.

1. Using your administrative login information, log into the administrative settings page.
2. Tap **System Portal Connection**.
3. On the next page, tap the **Upload logs** button.